


Inlet Supplies for Thermo TRACE 1300 GCs



Septa (11 mm)

Thermolite Septa	
Usable to 340 °C inlet temp.	
qty.	cat.#
50-pk.	27142
100-pk.	27143

Premium Non-Stick BTO Septa	
Usable to 400 °C inlet temp.	
qty.	cat.#
50-pk.	27090
100-pk.	27091

Viton® O-Rings

Description	Max Temp	10-pk. cat.#	50-pk. cat.#
Viton O-Rings	300 °C	22241	22242

Sky® Inlet Liner

ID x OD x L	qty.	cat.#
Single Taper, Sky Technology, Borosilicate Glass with Quartz Wool	ea.	23447.1
4.0 mm x 6.5 mm x 78.5 mm	5-pk.	23447.5

ID x OD x L	qty.	cat.#
Single Taper, Sky Technology, Borosilicate Glass	ea.	23446.1
4.0 mm x 6.5 mm x 78.5 mm	5-pk.	23446.5

ID x OD x L	qty.	cat.#
Straight, Sky Technology, Borosilicate Glass with Quartz Wool	ea.	23444.1
4.0 mm x 6.3 mm x 78.5 mm	5-pk.	23444.5

Dual Vespel® Ring Inlet Seals

0.8 mm ID Dual Vespel Ring Inlet Seal	2-pk. cat.#	10-pk. cat.#	50-pk. cat.#
Gold-Plated	21240	21241	23418

Replacement Inlet Seals With Washers

Single-Column 0.8 mm ID (Opening)	2-pk. cat.#	10-pk. cat.#	50-pk. cat.#
Gold-Plated	21317	21318	23415

RESTEK
Pure Chromatography